

BD GasPak™ EZ Gas Generating Systems

GasPak™ EZ Standard Incubation Container 15 Plates.	260671
GasPak™ EZ Large Incubation Container 30 Plates.....	260672
GasPak™ EZ Anaerobe Container System Sachets 20 Sachets Per Box.....	260678
GasPak™ EZ CO ₂ Container System Sachets 20 Sachets Per Box.....	260679
GasPak™ EZ Campy Container System Sachets 20 Sachets Per Box.....	260680
GasPak™ EZ Anaerobe Pouch System 20 Sachets, Resealable Pouches and Dry Anaerobic Indicators.....	260683
GasPak™ EZ CO ₂ Pouch System 20 Sachets and Resealable Pouches.....	260684
GasPak™ EZ Campy Pouch System 20 Sachets and Resealable Pouches.....	260685
GasPak™ EZ Standard Incubation Container Rack Holds 15-18 Petri Dishes.....	260673
GasPak™ EZ Large Incubation Container Rack Holds 30-33 Petri Dishes.....	260674
BBL™ GasPak™ Dry Anaerobic Indicator Strips Package of 100.....	271051
BBL™ GasPak™ Disposable Anaerobic Indicator Package of 100.....	270504
BBL™ GasPak™ 100 System.....	260626
BBL™ GasPak™ 150 Large Anaerobic System.....	260628

Complete Solutions for Anaerobic, Campylobacter and CO₂ Enriched Testing

- BD BBL™ Prepared Plated Media
 - CDC Anaerobe Agar
 - Brucella Agar
 - Campylobacter Agar with 5 Anti-microbics
- BD Vacutainer® Anaerobic Specimen Collector
- BD Port-A-Cul™ Anaerobic Transport Systems
- BBL™ CultureSwab™ Plus Collection and Transport Systems
- BBL™ QualiSwab™ Anaerobe Set
- BD GasPak™ EZ Gas Generating Systems
- BBL™ GasPak™ Environmental Generating Systems
- BD Bio-Bag™ Disposable Environmental Chambers
- Indicators
 - Anaerobic and CO₂ indicator strips



BD Diagnostics
Tullastrasse 8-12
69126 Heidelberg
Germany
Tel. (49) 6221 305 0
Fax (49) 6221 305 604
www.bd.com

BD, BD Logo and all trademarks are the property of Becton, Dickinson and Company. CultureSwab is a trademark of Becton, Dickinson and Company. ©2010 BD - XBA112540 Printed in Europe
Becton Dickinson GmbH - General Manager: Matthias Boos - Registered Office: Heidelberg - Commercial Register: Mannheim HRB 330 707

BD GasPak™ EZ Gas Generating System

Unbeatable Convenience in a Complete Integrated System



Helping all people
live healthy lives

A commitment
to innovation
and performance

In the late 1960s BD launched the first commercial system for environmental generation BD GasPak jars and envelopes. This innovative line eliminated the need for vacuum pumps, gas tanks, nanometers and measure release valves required at that time to create an appropriate growth environment for micro-organisms. Over the next 35 years BD continually introduced innovative products in response to our customers' evolving needs - including GasPak Plus envelopes, GasPak Pouch System, Bio-Bag™ Systems, anaerobic and CO₂ indicators and a complete line of accessories. Today, in keeping with our history of innovation, BD is pleased to present the GasPak EZ Gas Generating System. The GasPak EZ System represents the first and only complete, integrated system for rapid, waterless and catalyst-free gas generation.

BD GasPak™ EZ Gas Generating Systems

BD GasPak™ EZ Container



BD GasPak™ EZ Container Systems offer waterless, catalyst-free convenience for use in the innovative new GasPak EZ Incubation Containers or the original GasPak 100 or GasPak 150 jars. The GasPak EZ Gas Generating Container Sachets produce anaerobic, microaerophilic (reduced oxygen) or CO₂ enriched environments.


Easy as 1, 2, 3



1 Place inoculated BBL™ plates in the GasPak EZ Petri dish rack and lower rack into GasPak EZ Incubation Container.



2 Remove appropriate number of GasPak EZ Gas Generating Sachets and place in container. (Add BBL™ Anaerobic Indicator when using the GasPak EZ Anaerobic System.)



3 Latch GasPak EZ Incubation Container lid in place and incubate for culture.

Container System Features:

- Two GasPak EZ Incubation Container Sizes (15 and 30 Petri dish capacity)
 - > Containers are non-breakable and chemical resistant for maximum safety
 - > Easy-close latches assure a secure, airtight seal
 - > Removable Petri dish racks ease workflow
 - > Containers are designed to hold standard test tube racks for greater convenience
- GasPak EZ Sachets are endorsed for use in both the new GasPak EZ Containers and the original GasPak 100 and GasPak 150 jars
- Sachets require no water or catalyst for activation
- Sachets are engineered to provide rapid gas generation for optimal bacterial growth

BD GasPak™ EZ Pouch Systems

BD GasPak EZ Pouch Systems offer the convenience of pouches integrated into complete kits with everything you need to generate a pouch based anaerobic, microaerophilic or CO₂ enriched environment.

Easy as 1, 2, 3



1 Place inoculated plates inside resealable pouch.



2 Remove one GasPak EZ Gas Generating Sachet and place in pouch. (Add BBL™ Anaerobic Indicator when using the GasPak EZ Anaerobic System.)



3 Seal pouch by pressing zipper together and incubate for culture.

Pouch System Features:

- The only complete gas generation kit – pouches, sachets and indicators are all provided
- Large pouch capacity
 - > One to four Petri dish capacity for Anaerobe Pouch System
 - > One to two Petri dish capacity for CO₂ and Campy Pouch Systems
- One-step resealable pouch. No sealing bars needed!
- Custom multi-layer plastic pouch material preserves environmental integrity for maximum performance
- Sachets require no water or catalyst for activation
- Sachets are engineered to provide rapid gas generation

